

INSTRUCTIONS FOR USE of medical device

Injectable implant based on hyaluronic acid DIART

Composition

Injectable implant based on hyaluronic acid, 1.1 %, 1 ml, contains:

Sodium hyaluronate 11.0 mg
Succinate buffer pH 7.4..... up to 1.0 ml

Injectable implant based on hyaluronic acid, 1.1 %, 2 ml, contains:

Sodium hyaluronate 22.0 mg
Succinate buffer pH 7.4..... up to 2.0 ml

Injectable implant based on hyaluronic acid, 1.8 %, 1 ml, contains:

Sodium hyaluronate 18.0 mg
Succinate buffer pH 7.4..... up to 1.0 ml

Injectable implant based on hyaluronic acid, 1.8 %, 2 ml, contains:

Sodium hyaluronate 36.0 mg
Succinate buffer pH 7.4..... up to 2.0 ml

Contents of packaging

Injectable implant based on hyaluronic acid **DIART** in a pre-filled syringe with two needles.

Description

Hyaluronic acid is a natural polysaccharide (glycosaminoglycan) which is an important structural element of the connective tissue of the skin and synovial fluid.

DIART is a colorless, clear, elastic gel of non-stabilized hyaluronic acid of non-animal origin. Gel is sterile, pyrogen-free, with physiological pH.

Mode of action

DIART is a biologically compatible injectable long-acting implant with viscoelastic properties, acting as a buffer; the liquid that restores mobility of joints.

Indications

Temporary replacement of synovial fluid in synovial fluid deficiency in the knee joint, caused by injuries or degenerative diseases, such as osteoarthritis.


Contraindications

- Hypersensitivity to any component of medical device.
- Dysfunction of blood coagulation.
- Age under 18 years.

Method of use

DIART should be administered only by healthcare professionals with an experience in injecting the implants. Before the procedure, a physician should become familiar with the medical history of the patient and inform her/him about the projected therapeutic outcomes and possible undesirable effects.

The physician should attach the label with the product data (batch number and expiry date) to the Medical Record of the patient undergoing the procedure.



During intra-articular injections, the physician should strictly follow the aseptic technique. It is recommended to inject the implant using the needles that are in the pack. If the needle passage is impaired and injection of the implant is complicated, stop the procedure and replace the needle. If applicable, intra-articular injection is recommended to be US-guided. **DIART** should be injected into the affected joint at the dose of 1–2 ml depending on the size of the joint. Treatment course consists of 3 injections into one joint at 1-week intervals. If required, a repeated course may be administered.

Adverse reactions

Hypersensitivity reactions, pulsing sensation in the joints, myalgia. Local reactions: flushing, edema, bruising, itching, minor pain in the injection site. Individual reactions caused by the invasive procedure are possible such as dizziness, nausea, headache, gastric disturbance, limited joint movement. These undesirable effects are temporary and usually resolve within a few hours after the injection.

Joint infections and septic arthritis may occur in case of a non-compliance with aseptic conditions during the procedure.


Due to the injection of the medical device in the presence of inflammatory process and inflammation of the synovial membrane of the joint, there may be increased pain, the development of aseptic inflammation or exacerbation of symptoms of exudative inflammation of the synovial membrane.

When the medical device is injected into the meniscus, fat body, intra-articular ligaments, acute pain or significant exacerbation of the pain syndrome may be observed.

In the event of any adverse reactions, the use of the medical device should be discontinued immediately and the physician and the manufacturer should be informed.

Limitations, precautions and warnings

- No clinical data is available regarding the tolerance of **DIART** by pregnant or breastfeeding women, thus **DIART** should not be used in women who are pregnant or breastfeeding.
- Do not use with other routes of administration, except for intra-articular.
- Within a week before the procedure, the patients should not take medicines that may affect blood clotting function (e.g., acetylsalicylic acid, anticoagulants). This may cause hemorrhagic complications during the procedure. When necessary, for example, if the patient is treated with medicines of such type, a doctor should be consulted.
- Injection of the medical device in the areas of the nerves and their plexuses requires extreme caution. Impaired technique of implant injection in these areas may result in traumatic injury to the peripheral nerves, increased pain syndrome, sensory impairment, and the development of paralysis in the appropriate regions of innervation.
- Avoid injection of the medical device into the tenosynovial membranes of the peroneal, tibial or other adjacent muscles, which may lead to increased pain syndrome, motor disorders of the damaged joint and adjacent joints of the leg.
- It is recommended to avoid any excessive physical exercise on the joints within 48 hours after procedure (such as playing tennis, jumping, football, etc.).
- The injection site should not be exposed to intense heating or cooling during the first 48 hours after injection.
- To prevent skin itching, the use of alcohol-free disinfectants is recommended.
- The syringe, needles and the residues of unused medical device are to be disposed immediately after the procedure.
- Do not mix with other devices.
- Prior to application, verify the package integrity and check its shelf life. Do not use the device after the expiration date or if the package integrity is damaged.
- For single use only. Repeated use can lead to infection. Do not re-use.



– **DIART** is NOT the subject to resterilization.

Note: this medical device should be only used by the experienced staff at the designated clinics and in accordance with the current regulations.

Incompatibilities

Injectable implants **DIART** should not be co-administered with quaternary ammonium salts (benzalkonium chloride) and chlorhexidine.

Storage conditions

Store in a dry, dark place at +5 °C to +30 °C. Keep away from children. Do not freeze.

Shelf life: 2 years. Shelf life is valid if the storage conditions are observed and the package is not damaged.

Authorized representative in the European Community

Diaco Biofarmaceutici S.R.L.
Via Flavia, 124, 34147, Trieste, Italy.
E-mail: info@diaco.it
<https://www.diaco.it>


Name and address of manufacturer

Yuria-pharm LLC, 10, M. Amosova Str., Kyiv, Ukraine, 03038.
Tel.: +38 (044) 275-92-42, +38 (044) 275-01-08.
E-mail: uf@uf.ua
<https://www.uf.ua>
Manufacturing site address: 108, Kobzarska Str., Cherkasy, Ukraine, 18030.

















If you have any comments on the medical device or would like to give us feedback, please use the following options to contact us:

- 1) email us at feedback@uf.ua;
- 2) send a text message via Viber, Telegram or WhatsApp to the number: +38 (095) 275-33-01;
- 3) call us at +38 (095) 275-33-01 or +38 (0800) 401-771 (charged in accordance with your operator's tariff plan).



Graphical symbols and their interpretation	
	Sterilized using steam or dry heat
	Sterile medical device in primary packaging
	Do not re-use
	Caution
	Consult instructions for use

	Manufacturer
	Keep away from sunlight
	Keep dry
	Do not resterilize
	Do not use if package is damaged
	Fragile, handle with care
	Non-pyrogenic
	Temperature limit
	Date of manufacture
	Use-by date
	Batch code
 UA.TR.099	Conformity mark to the technical regulations and the code of the conformity assessment body/Знак відповідності технічним регламентам та код органу з оцінки відповідності
	Mark of compliance with Directive 93/42/EEC on medical devices and the Notified Body
	Authorized representative in the European Community

Date of last revision of instructions for use: 07.03.2023.
Version: 07.

